Maria Cristina Mansur, M.D.

GentleLASE® ile Çillerin Tedavisi

Dr. Maria Cristina Mansur,
Tıbbi Direktör, Dermatoloji
Bölümü Centro Científico de
Medicina Cutânea—MG;
Professor, Universidade Federal de
Juiz de Fora—MG, Brezilya

Giriş

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Referanslar

1. Rosenbach A; Lee SJ, Johr RH. "Treatment of medium-brown solar lentigines using an alexandrite laser designed for hair reduction." Arch Dermatol 138 (4):547-8, 2002.

2. Rashid T, Hussain I, Haider M, Harron TS. "Laser therapy of freckles and lentigines with quasicontinuous, frequency-doubled, Nd:YAG (532 nm) laser in Fitzpatrick skin type IV: A 24 month follow up." J Cosmet Laser Ther 2002; 4:81-85.

3. Kawada A, Shiraishi H, Asai M et al. "Clinical improvement of solar lentigines and ephelides with an intense pulsed light source." Dermatol Surg 2002; 28:504-8.

4. Chan HH, Murad A, Kono T, Dover J. "Clinical application of lasers in Asians." Dermatol Surg 2002; 28:556-63

5. Lam AY, Wong DS, Lam LK, Ho WS, Chan HH. "A retrospective study on the efficacy and complications of Q-switched alexandrite laser in the treatment of acquired bilateral nevus of Ota-like macules." Dermatol Surg 2001; 27:937-41.

6. Somoyos K, Boonchu K, Somsak K, Pannadda L, Leopairut J. "Copper vapour laser treatment of cafe-au-lait macules." Br J Dermatol 1996; 135:964-96.

7. Chan HH, Fung WKK, Ying SY, et al. "An in vivo trial comparing the use of different types of 532 nm Neodymium: Yttrium-Aluminum- Garnet (Nd:YAG) lasers in the treatment of facial lentigines in oriental patients." Dermatol Surg 2000; 26:743-749

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Candela Corporation
530 Boston Post Road
Wayland, MA 01778, USA
Phone: (508) 358-7637
Fax: (508) 358-5569
Toll Free: (800) 821-2013
www.candelalaser.com

CANDELA
BE CONFIDENT